


### SIZE INFORMATION

EzyAs™ applicators are available in three sizes and colour coded for easy identification.



## APPLICATOR SIZE GUIDE MEASUREMENT A OR B

CENTIMETRES	INCHES	RECOMMENDED SIZE
28-35	11-13.5	SMALL
32-40	12.5-15.75	MEDIUM
38+	15+	LARGE

Available from:



International: +61 418 386 637 Email: info@ezyasabc.com www.ezyasabc.com Episode & People's
Choice Winner'
Plus
Grand Final
Les is More - Big Heart
Award Winner'
the New Inventors

Failure to adhere to care instructions will void the lifetime guarantee



# GARMENT APPLICATOR FOR COMPRESSION STOCKINGS





EzyAs<sup>™</sup> is a unique device designed to assist in the application of compression stockings and garments.

EzyAs<sup>™</sup> enables a wide range of compression materials to be applied over a patient's leg or arm without drag against fragile skin or disturbing dressings.



### **BENEFITS**

- © QUICK AND EASY TO USE
  WITH COMFORT AND CONTROL
- © PRACTICAL FOR EITHER LEG OR ARM APPLICATIONS
- MINIMAL DISRUPTION TO
  AN INJURY OR DRESSING
- COMPACT, LIGHTWEIGHT
  AND TRAVEL FRIENDLY
- THREE COLOUR CODED SIZES FOR EASY IDENTIFICATION
- © WORKS WITH OPEN AND CLOSED TOE STOCKINGS



### **TESTIMONIALS**

Health professionals confirm EzyAs™ applicators provide beneficial outcomes.

"Ezy-As reduces work load and the difficulty of donning stockings for the nursing staff."

Occupational Therapist, ACT

"Patients feel the Ezy-As applicator does not traumatise their skin integrity."

Woundcare Consultant, QLD

### HANDLE ATTACHMENT

The optional handle attachment\* has allowed many customers who have limited flexibility to be independent in applying their compression stockings.

The handle attachment\* is designed to fit any size applicator in the range.

\*The handle attachment is sold separately.



Visit us at www.ezyasabc.com to view instructional videos